

SIAMDENT.
Morita Group Company




SX3000

SIAMDENT.

Siamdent Co., Ltd.
Factory : 71/10 Mu 5, Thakham, Bangpakong, Chachuengsao 24130 Tel : 66 38 573042, 573133 Fax : 66 38 573043
Office : 444 Olympia Thai Tower, 3rd Floor, Ratchandapisek Road, Samsennok, Huay Kwang, Bangkok 10310 Tel. 02-512-6049-50 Fax. 02-512-6099



Standard features & Function




Assistant's Syringe Holder
High power suction syringe, saliva ejector and 3-way syringe are provided with multi-functional keypad switch



Spittoon Unit Body
Hygienic ceramic spittoon with cup filler, water basin and assistant's instrument holders




LED Operating Light
Bright shadow-free illumination & detachable silicon handle cover.



Extendable
Adjustable positioning for children & adults



Water Tank
Distilled water is supplied for handpiece



Doctor Tray Table
Spacious tray table with 5 instrument holders, multi-functional keypad switch dental size film viewer and First priority system



Armrest
Well-balanced single armrest provided



Base Safety
Emergency stop under base cover



Junction Box
Convenient for installing the chair unit



Foot Controller
Easy-to-use one foot pedal with multi-functions

Standard Equipment & Feature

Doctor side

- * Doctor tray table with 5 instrument holders
- * High speed air turbine circuit with fiber optic 4H tubing (2)
- * Low speed airmotor circuit with 4H tubing (1)
- * Doctor side 3-way syringe with tubing (1)
- * Doctor side multi-functional keypad switch
- * Dental size film viewer

Assistant Side

- * Assistant instrument 3 holders
- * High power suction syringe
- * Saliva ejector on assistant side
- * 3-way syringe with tubing
- * Assistant side multi-functional keypad switch

Spittoon Unit Body

- * Hygienic ceramic spittoon with water basin
- * Sensor-type cup filler
- * Drain cut-filter
- * Water tank
- * Manual water control knobs for water basin
- * On-Off switches for basin water & cup filler

LED Operating light

- * On-Off by automatic
- * Operating light On-Off by sensor
- * Dimmer light

Foot Controller

- * Foot pedal for handpieces & levers for chair movements
- * Switch for water supply for handpieces
- * Switch for light On-Off
- * One preset positioning & zero position
- * Emergency stop by pedal

Armrest

- * Single armrest

Junction box

- * Main electric, air & water switches
- * Adjusted to local line voltage

Chair Driving System

- * Hydraulic system
- * Two preset positioning & zero position
- * Gargling positioning & last treatment positioning

Vacuum System

- * Air vacuum system or motor vacuum system

Optional Equipment & Features

- * Micromotor circuit & tubing
- * 2H, 3H tubings for high & low speed
- * Built-in type scaler
- * Panorama & Ceph size film viewer
- * LCD monitor support
- * Intraoral camera

- * Twin armrest instead of single armrest
- * Motor suction vacuum system wet or dry type
- * Water venturi saliva ejector
- * Water Warmer
- * Emergency stop under base cover
- * Monitor support

We reserve the right to change the specifications without prior notice

Collor Guide



Pink



Orange



Green



Blue



Violet